

**B-172 MYELOPEROXIDASE STAIN USING  
4-CHLORO-1-NAPHTHOL**

**FIXATION:** Streak thin (approx. one cell thick) smears from fresh non-anticoagulated blood across a sterile slide by means of a second slide or cover glass. Fix in 10% Formal-Alcohol for 60 seconds.

**PREPARATION OF WORKING SOLUTIONS:**

A. 4-Chloro-1-Naphthol Working Solution

Just before use mix:

- |   |      |
|---|------|
| 1. <b><u>4-Chloro-1-Naphthol Stock Solution</u></b> (B-172-1) | 2ml  |
| 2. <b><u>Tris Buffer, pH 7.4</u></b> (B-172-2)                | 38ml |
| 3. <b><u>Hydrogen Peroxide, 0.15%</u></b> (see below)         | 1ml  |

B. Hydrogen Peroxide, 0.15%

Dilute 1ml of **Hydrogen Peroxide, 3%** (B-172-3A), to 20ml (qs.) with distilled water and mix well. Prepare fresh before use.

**STAINING PROCEDURE:**

1. Wash slides with gently running tap water for 15 seconds and shake off excess.
2. Incubate the smears in **4-Chloro-1-Naphthol Working Solution** (see directions above) for 8-10 minutes at room temperature.
3. Wash the smears with distilled water.
4. **OPTIONAL.** Counterstain for 1 to 3 minutes with **Methyl Green, 1% Aqueous** (B-172-4) or use a **Wright-Giemsa** or **Giemsa Method** (such as B-153).
5. Rinse the smears rapidly with cold distilled water and blot dry.

**RESULTS:**

Myeloblast and Neutrophil Cytoplasm.....blue-black granules

The negative unstained image of the nucleus is easily discerned, monocytes are moderately stained but all other formed elements of the blood in both bone-marrow and peripheral blood films do not stain for myeloperoxidase.

**REFERENCES:**

Elias, J. M., American Journal of Clinical Pathology, Vol. 73, No. 6, June 1980  
Omri M. B., State University of New York, 3.0-006