

F-400 AYOUB-SHKLAR METHOD FOR KERATIN & PREKERATIN

FIXATION: Buffered Neutral Formalin, 10% (F—113)

SECTIONS: Paraffin @ 6 microns

STAINING:

1. Deparaffinize and hydrate to distilled water.
2. Acid Fuchsin Solution (F-400-1) for 3 minutes.
3. Transfer slides directly to Aniline Blue-Orange G Solution (F-400-2) for 45 minutes.
4. Transfer slides directly to 95% alcohol for three changes.
5. Dehydrate with absolute alcohol and clear in Xylene (C-120), two changes each.
6. Mount with Permout (M-18).

STAINING RESULTS:

Keratin	Brilliant Red
Stratified squamous epithelium	Gray
Stratum spinosum, upper area	
Illustrating prekeratin-like	
Substances	Orange
Connective Tissue	Deep Blue
Erythrocytes	Deep Red

REFERENCES:

AFIP Manual of Histological Staining Methods, 3rd ed., Ed. L. Luna:
New York: McGraw Hill Publications, c. 1968, p. 82.

Ayoub, P. and Shklar, G., J. Oral Surg., 16:580-581, 1963.