




CD117, c-kit

USA IVD

This is able to differentiate gastrointestinal stromal tumors from Kaposi's sarcoma and tumors of smooth muscle origin.

Rabbit Monoclonal. Clone: YR145

Control: GIST. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Cytoplasmic, Membranous



GIST

Catalog no.	Volume
C17-2010-1	0.1 ml concentrated
C17-2010-2	0.5 ml concentrated
C17-2010-3	1ml concentrated
C17-2010-4	1ml Prediluted
C17-2010-5	7 ml Prediluted
C17-2010	5 Positive control slides