

## **L-757, WEIGERT'S IRON HEMATOXYLIN**

**Fixation:** Any general Fixative

**Sections:** Paraffin

**Staining:**

1. Deparaffinize and hydrate to distilled water.
2. Stain in working **Weigert's Hematoxylin** [equal volumes of **Part A** (L757-1) and **Part B** (L-757-2) for 5-6 minutes]
3. Counterstain as desired. **Van Gieson's** (F-757-3) for 3-5 minutes is recommended. Rinse in distilled water.
4. Differentiate in 95% alcohol. Dehydrate, clear and mount.

**Staining Results:**

Nuclei	Deep blue to black
--------	--------------------

**References:**

Conn, H.J.: Biological Stains, Williams & Wilkins Co., Balt., 1969, 8<sup>th</sup> ed., p.436.

Clark, G.: Staining Procedures, Williams & Wilkins Co., Balt., 1973, 3<sup>rd</sup>. ed., p. 38