




A-1-Antitrypsin

USA IVD

May be used in the detection of benign and malignant lesions of a histiocytic nature.

Rabbit Polyclonal

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Tonsil

Catalog no.	0.1 ml concentrated
AA-2010-1	0.5 ml concentrated
AA-2010-2	1 ml concentrated
AA-2010-3	1 ml Prediluted
AA-2010-4	7 ml Prediluted
AA-2010-5	5 Positive control slides
AA-2010	

Page

2

Volume