




ACTH

USA IVD

This is a useful marker in classification of pituitary tumors and the study of pituitary disease.

Rabbit Polyclonal

Control: Pituitary. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Pituitary

Catalog no.	Volume
AC-2010-1	0.1 ml concentrated
AC-2010-2	0.5 ml concentrated
AC-2010-3	1 ml concentrated
AC-2010-4	1 ml Prediluted
AC-2010-5	7 ml Prediluted
AC-2010	5 Positive control slides