




# Actin, Smooth Muscle

USA IV D

This will stain myofibroblasts and myoepithelial cells.

**Mouse Monoclonal, Clone: 1A4**

**Control:** Appendix. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG/K **Cellular localization:** Cytoplasmic



Appendix

<b>Catalog no.</b>	0.1 ml concentrated
ASM-2010-1	0.5 ml concentrated
ASM-2010-2	1 ml concentrated
ASM-2010-3	1 ml Prediluted
ASM-2010-4	7 ml Prediluted
ASM-2010-5	5 Positive control slides
ASM-2010	

**Volume**