




Alpha-Fetoprotein

USA IVD

Positive staining with this is seen in hepatocytes of fetal liver and hepatoma.

Rabbit Polyclonal

Control: Fetal liver. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Cellular localization:** Cytoplasmic



Fetal Liver

Catalog no.	Volume
AF-2010-1	0.1 ml concentrated
AF-2010-2	0.5 ml concentrated
AF-2010-3	1 ml concentrated
AF-2010-4	1 ml Prediluted
AF-2010-5	7 ml Prediluted
AF-2010	5 Positive control slides