




# Androgen Receptor

USA IVD

This has been useful in differentiating morpheaform basal cell carcinoma from desmoplastic trichoepithelioma in the skin.

**Mouse Monoclonal, Clone: AR441**

**Control:** Prostate. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG1 **Cellular localization:** Nuclear



Prostate

Catalog no.	Volume
AR-2010-1	0.1 ml concentrated
AR-2010-2	0.5 ml concentrated
AR-2010-3	1 ml concentrated
AR-2010-4	1 ml Prediluted
AR-2010-5	7 ml Prediluted
AR-2010	5 Positive control slides