



# Annexin A1


USA IVD

This is strongly expressed on the cell membrane and occasionally in the cytoplasm of tumor cells in 97% of samples from patients with hairy cell leukemia.

**Mouse Monoclonal, Clone: MRQ-3**

**Control:** Hairy Cell Leukemia. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

**Shelf life:** Up to 36 months. **Isotype:** IgG1 **Cellular localization:** Cytoplasmic, Membranous



Hairy Cell Leukemia

<b>Catalog no.</b>	<b>Volume</b>
ANN-2010-1	0.1 ml concentrated
ANN-2010-2	0.5 ml concentrated
ANN-2010-3	1 ml concentrated
ANN-2010-4	1 ml Prediluted
ANN-2010-5	7 ml Prediluted
ANN-2010	5 Positive control slides