



## BG8,Lewis


USA IV D

This has the ability to distinguish between pulmonary adenocarcinoma and epithelioid mesothelioma.

**Mouse Monoclonal, Clone: F3**

**Control:** Adenocarcinoma of lung. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

**Shelf life:** Up to 36 months. **Isotype:** IgM **Cellular localization:** Cytoplasmic



Adenocarcinoma of lung

### Catalog no.

### Volume

BL-210-1	
BL-2010-2	0.1 ml concentrated
BL-2010-3	0.5 ml concentrated
BL-2010-4	1 ml concentrated
BL-2010-5	1 ml Prediluted
BL-2010	7 ml Prediluted
	5 Positive control slides

