




BOB.1

USA IVD

This is a critical determinant of octamer-dependent gene transcription in B-lymphocytes.

Mouse Monoclonal, Clone: MRQ-35

Control: B-cell I lymphoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgM **Cellular localization:** Cytoplasmic, Nuclear



B-cell lymphoma

Catalog no.	Volume
BO-210-1	
BO-2010-2	0.1 ml concentrated
BO-2010-3	0.5 ml concentrated
BO-2010-4	1 ml concentrated
BO-2010-5	1 ml Prediluted
BO-2010	7 ml Prediluted
	5 Positive control slides