




## C4d

USA IVD

This has been a significant predictor of transplant kidney graft survival and as an aid in treating acute rejection.

### Rabbit Polyclonal

**Control:** Lymph node, Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.  
**Shelf life:** Up to 36 months. **Cellular localization:** Membranous, Cytoplasmic.



Lymph node

Catalog no.	Volume
C4-2010-1	0.1 ml concentrated
C4-2010-2	0.5 ml concentrated
C4-2010-3	1 ml concentrated
C4-2010-4	1 ml Prediluted
C4-2010-5	7 ml Prediluted
C4-2010	5 Positive control slides