



## CA-125


USA IVD

This recognizes specific carbohydrate-associated epitopes localized in the variable region of the cancer cell-expressed immunoglobulin heavy chains.

**Mouse Monoclonal, Clone: OC125**

**Control:** Ovarian Carcinoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

**Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Cytoplasmic, Membranous



Ovarian Carcinoma

### Catalog no.

### Volume

C12-2010-1	0.1 ml concentrated
C12-2010-2	0.5 ml concentrated
C12-2010-3	1 ml concentrated
C12-2010-4	1 ml Prediluted
C12-2010-5	7 ml Prediluted
C12-2010	5 Positive control slides