




CA15-3

USA IVD

This has been used for evaluating the primary site of a metastatic carcinoma of unknown origin and distinguishing between benign and malignant lesions.

Mouse Monoclonal, Clone: DF3

Control: Breast. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic



Breast

Catalog no.	Volume
C15-2010-1	
C15-2010-2	0.1 ml concentrated
C15-2010-3	0.5 ml concentrated
C15-2010-4	1 ml concentrated
C15-2010-5	1 ml Prediluted
C15-2010	7 ml Prediluted
	5 Positive control slides