




CD10

USA IVD

This is a useful marker for the characterization of childhood leukemia and B-cell lymphomas.

Mouse Monoclonal, Clone: 56C6

Control: Tonsil, Kidney. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Membranous, Cytoplasmic



Tonsil

Catalog no.	0.1 ml concentrated
C10-2010-1	0.5 ml concentrated
C10-2010-2	1 ml concentrated
C10-2010-3	1 ml Prediluted
C10-2010-4	7 ml Prediluted
C10-2010-5	5 Positive control slides
C10-2010	