




CD138

USA IV D

This is a marker for labeling normal and neoplastic plasma cells as well as plasmacytoid lymphoma, multiple myeloma, lymphoplasmacytic lymphoma, and various forms of Hodgkin's disease.

Mouse Monoclonal, Clone: B-A38

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ Cellular **localization:** Cytoplasmic, Membranous



Tonsil

Catalog no.

Volume

| | |
|------------|---------------------------|
| C38-2010-1 | 0.1 ml concentrated |
| C38-2010-2 | 0.5 ml concentrated |
| C38-2010-3 | 1ml concentrated |
| C38-2010-4 | 1ml Prediluted |
| C38-2010-5 | 7 ml Prediluted |
| C38-2010 | 5 Positive control slides |