




CD138

USA IV D

This is a marker for labeling normal and neoplastic plasma cells as well as plasmacytoid lymphoma, multiple myeloma, lymphoplasmacytic lymphoma, and various forms of Hodgkin's disease.

Mouse Monoclonal, Clone: B-A38

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ Cellular **localization:** Cytoplasmic, Membranous



Tonsil

Catalog no.	Volume
C38-2010-1	0.1 ml concentrated
C38-2010-2	0.5 ml concentrated
C38-2010-3	1ml concentrated
C38-2010-4	1ml Prediluted
C38-2010-5	7 ml Prediluted
C38-2010	5 Positive control slides