




# CD163

USA IVD

This is helpful in distinguishing synovial macrophages from synovial intimal fibroblasts in the setting of rheumatoid arthritis.

**Mouse Monoclonal, Clone: MRQ-26**

**Control:** Inflamed Tissue. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic, Membranous



Inflamed Tissue

<b>Catalog no.</b>	0.1 ml concentrated
C16-2010-1	0.5 ml concentrated
C16-2010-2	1ml concentrated
C16-2010-3	1ml Prediluted
C16-2010-4	7 ml Prediluted
C16-2010-5	5 Positive control slides
C16-2010	

## Volume