




CD19

USA IVD

This is a reliable surface marker of this lineage over a wide range of maturational stages.

Mouse Monoclonal, Clone: MRQ-36

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Membranous



Tonsil

Catalog no.	Volume
C19-2010-1	0.1 ml concentrated
C19-2010-2	0.5 ml concentrated
C19-2010-3	1 ml concentrated
C19-2010-4	1 ml Prediluted
C19-2010-5	7 ml Prediluted
C19-2010	5 Positive control slides