




## CD1a

USA IV D

This has been used to differentiate various cutaneous lymphomas (T-cell) from B-cell lymphomas and pseudolymphomas.

**Mouse Monoclonal, Clone: MTB1**

**Control:** Thymus. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8<sup>0</sup>C. **Shelf life:** Up to 36 months. **Isotype:** IgG1/K **Cellular localization:** Membranous



Thymus

Catalog no.	Volume
CD1-2010-1	
CD1-2010-2	0.1 ml concentrated
CD1-2010-3	0.5 ml concentrated
CD1-2010-4	1 ml concentrated
CD1-2010-5	1 ml Prediluted
CD1-2010	7 ml Prediluted
	5 Positive control slides

