




## CD21

USA IVD

This is useful in the identification of follicular dendritic cell matrix found in normal lymph nodes and tonsillar tissue.

**Rabbit Monoclonal, Clone: EP3093**

**Control:** Tonsil . **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** 2G9, IgG<sub>2a</sub>, EP3093, IgG **Cellular localization:** Membranous, Cytoplasmic



Tonsil

### Catalog no.

	<b>Volume</b>
C21-2010-1	
C21-2010-2	0.1 ml concentrated
C21-2010-3	0.5 ml concentrated
C21-2010-4	1 ml concentrated
C21-2010-5	1 ml Prediluted
C21-2010	7 ml Prediluted
	5 Positive control slides

