



## CD23

USA IVD

This is used to differentiate between B-CLL and B-SII's from mantle cell lymphomas and small cleaved lymphomas.

**Mouse Monoclonal, Clone: 1B12**

**Control:** Tonsil . **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Membranous



Tonsil

Catalog no.	Volume
C23-2010-1	0.1 ml concentrated
C23-2010-2	0.5 ml concentrated
C23-2010-3	1 ml concentrated
C23-2010-4	1 ml Prediluted
C23-2010-5	7 ml Prediluted
C23-2010	5 Positive control slides