




CD25

USA IVD

This is used to distinguish neoplastic mast cell aggregates from reactive proliferations.

Mouse Monoclonal, Clone: 4C9

Control: Lesions of Mastocytosis. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2b} **Cellular localization:** Membranous, Cytoplasmic



Lesions of Mastocytosis

Catalog no.	0.1 ml concentrated
C25-2010-1	0.5 ml concentrated
C25-2010-2	1 ml concentrated
C25-2010-3	1 ml Prediluted
C25-2010-4	7 ml Prediluted
C25-2010-5	5 Positive control slides
C25-2010	

Volume