




## CD25

USA IVD

This is used to distinguish neoplastic mast cell aggregates from reactive proliferations.

**Mouse Monoclonal, Clone: 4C9**

**Control:** Lesions of Mastocytosis. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>2b</sub> **Cellular localization:** Membranous, Cytoplasmic



Lesions of Mastocytosis

<b>Catalog no.</b>	0.1 ml concentrated
C25-2010-1	0.5 ml concentrated
C25-2010-2	1 ml concentrated
C25-2010-3	1 ml Prediluted
C25-2010-4	7 ml Prediluted
C25-2010-5	5 Positive control slides
C25-2010	

**Volume**