




CD3

USA IVD

This is the best all around T-cell marker.

Rabbit Polyclonal

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Cellular localization:** Membranous



Tonsil

Catalog no.	Volume
CD3-2010-1	
CD3-2010-2	0.1 ml concentrated
CD3-2010-3	0.5 ml concentrated
CD3-2010-4	1 ml concentrated
CD3-2010-5	1 ml Prediluted
CD3-2010	7 ml Prediluted
	5 Positive control slides