



CD30


USA IVD

This detects a formalin-resistant epitope that is expressed by Reed Sternberg cells.

Mouse Monoclonal, Clone: Ber-H2

Control: Hodgkin's Lymphoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C.

Shelf life: Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Membranous



Hodgkin's Lymphoma

Catalog no.	Volume
C30-2010-1	0.1 ml concentrated
C30-2010-2	0.5 ml concentrated
C30-2010-3	1 ml concentrated
C30-2010-4	1 ml Prediluted
C30-2010-5	7 ml Prediluted
C30-2010	5 Positive control slides