




CD31

USA IVD

This is used to identify and concentrate cells of the hematopoietic system for experimental studies and bone marrow transplantation.

Mouse Monoclonal, Clone: JC70

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** JC70, IgG₁/K, 1A10, IgG₁/K. **Cellular localization:** Cytoplasmic, Membranous.



Tonsil

Catalog no.	Volume
C31-2010-1	0.1 ml concentrated
C31-2010-2	0.5 ml concentrated
C31-2010-3	1 ml concentrated
C31-2010-4	1 ml Prediluted
C31-2010-5	7 ml Prediluted
C31-2010	5 Positive control slides

