




CD34

USA IV D

This recognizes a cell surface antigen of 110kD that is expressed on human hematopoietic progenitor cells.

Mouse Monoclonal, Clone: QBEnd/10

Control: Tonsil, Placenta. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Membranous



Tonsil

Catalog no.	Volume
C34-2010-1	0.1 ml concentrated
C34-2010-2	0.5 ml concentrated
C34-2010-3	1 ml concentrated
C34-2010-4	1 ml Prediluted
C34-2010-5	7 ml Prediluted
C34-2010	5 Positive control slides