




CD35

USA IVD

This is a B-cell marker which labels follicular dendritic reticulum cells and tumors derived from such cells as follicular dendritic cell tumors.

Mouse Monoclonal, Clone: RLB25

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2b} **Cellular localization:** Membranous



Tonsil

Catalog no.	Volume
C35-2010-1	0.1 ml concentrated
C35-2010-2	0.5 ml concentrated
C35-2010-3	1 ml concentrated
C35-2010-4	1 ml Prediluted
C35-2010-5	7 ml Prediluted
C35-2010	5 Positive control slides