




CD4

USA IVD

This is used in the immunophenotyping of reactive lymphocytes and lymphoproliferative disorders.

Rabbit Monoclonal, Clone: SP35

Control: Tonsil . **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Membranous, Cytoplasmic



Tonsil

| Catalog no. | Volume |
|-------------|---------------------------|
| CD4-2010-1 | 0.1 ml concentrated |
| CD4-2010-2 | 0.5 ml concentrated |
| CD4-2010-3 | 1 ml concentrated |
| CD4-2010-4 | 1 ml Prediluted |
| CD4-2010-5 | 7 ml Prediluted |
| CD4-2010 | 5 Positive control slides |

