




CD44

USA IVD

This can differentiate between urothelial transitional cell carcinoma-in-situ from non-neoplastic changes in the urothelium.

Mouse Monoclonal, Clone: MRQ-13

Control: Benign Urothelium. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2a} Cellular **localization:** Membranous



Benign Urothelium

Catalog no.	0.1 ml concentrated
C44-2010-1	0.5 ml concentrated
C44-2010-2	1 ml concentrated
C44-2010-3	1 ml Prediluted
C44-2010-4	7 ml Prediluted
C44-2010-5	5 Positive control slides
C44-2010	

Volume