




CD45 (LCA)

USA IVD

This is used to aid the differential diagnosis of undifferentiated neoplasms, whenever malignant lymphoma is suspected by the morphological or clinical data.

Mouse Monoclonal, Clone: 2B11 & PD7/26

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Membranous



Tonsil

Catalog no.	0.1 ml concentrated
LCA-2010-1	0.5 ml concentrated
LCA-2010-2	1ml concentrated
LCA-2010-3	1ml Prediluted
LCA-2010-4	7 ml Prediluted
LCA-2010-5	5 Positive control slides
LCA-2010	