




## CD45R

USA IVD

This acts as a phosphor-tyrosine phosphatase, a vital component for the efficient tyrosine phosphorylation induction by TCR/CD3 complex.

**Mouse Monoclonal, Clone: MB1**

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> Cellular **localization:** Membranous



Tonsil

<b>Catalog no.</b>	0.1 ml concentrated
45R-2010-1	0.5 ml concentrated
45R-2010-2	1ml concentrated
45R-2010-3	1ml Prediluted
45R-2010-4	7 ml Prediluted
45R-2010-5	5 Positive control slides
45R-2010	