




CD45R

USA IVD

This acts as a phosphor-tyrosine phosphatase, a vital component for the efficient tyrosine phosphorylation induction by TCR/CD3 complex.

Mouse Monoclonal, Clone: MB1

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ Cellular **localization:** Membranous



Tonsil

Catalog no.	0.1 ml concentrated
45R-2010-1	0.5 ml concentrated
45R-2010-2	1ml concentrated
45R-2010-3	1ml Prediluted
45R-2010-4	7 ml Prediluted
45R-2010-5	5 Positive control slides
45R-2010	

Volume