




CD5

USA IVD

This is a T-cell marker that also reacts with a range of neoplastic B-cells.

Rabbit Monoclonal, Clone: SP19

Control: Tonsil . **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** SP19, IgG, 4V7, IgG/K **Cellular localization:** Membranous



Tonsil

Catalog no.	Volume
CD5-2010-1	0.1 ml concentrated
CD5-2010-2	0.5 ml concentrated
CD5-2010-3	1 ml concentrated
CD5-2010-4	1 ml Prediluted
CD5-2010-5	7 ml Prediluted
CD5-2010	5 Positive control slides