




CD56

USA IVD

This recognizes two proteins of the neural cell adhesion molecule.

Mouse Monoclonal, Clone: 123C3.D5

Control: Neuroblastoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K Cellular **localization:** Membranous



Neuroblastoma

Catalog no.	Volume
C56-2010-1	0.1 ml concentrated
C56-2010-2	0.5 ml concentrated
C56-2010-3	1ml concentrated
C56-2010-4	1ml Prediluted
C56-2010-5	7 ml Prediluted
C56-2010	5 Positive control slides