




## CD57

USA IVD

This marks a subset of lymphocytes, natural killer cells.

**Mouse Monoclonal, Clone: NK-1**

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgM/K Cellular **localization:** Membranous



Tonsil

<b>Catalog no.</b>	0.1 ml concentrated
C57-2010-1	0.5 ml concentrated
C57-2010-2	1ml concentrated
C57-2010-3	1ml Prediluted
C57-2010-4	7 ml Prediluted
C57-2010-5	5 Positive control slides
C57-2010	