




CD61

USA IVD

This is useful in identifying megakaryoblastic differentiation as seen in megakaryoblastic leukemia.

Mouse Monoclonal, Clone: 2f2

Control: Bone Marrow. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG1 **Cellular localization:** Cytoplasmic



Bone Marrow

Catalog no.	Volume
C61-2010-1	0.1 ml concentrated
C61-2010-2	0.5 ml concentrated
C61-2010-3	1ml concentrated
C61-2010-4	1ml Prediluted
C61-2010-5	7 ml Prediluted
C61-2010	5 Positive control slides