




CD63

USA IVD

This is useful in identifying malignant melanoma

Mouse Monoclonal, Clone: NK1/C3

Control: Melanoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic, Membranous



Melanoma

Catalog no.	Volume
C63-2010-1	0.1 ml concentrated
C63-2010-2	0.5 ml concentrated
C63-2010-3	1ml concentrated
C63-2010-4	1ml Prediluted
C63-2010-5	7 ml Prediluted
C63-2010	5 Positive control slides