




CD68

USA IVD

This can be useful for the identification of myelomonocytic and histiocytic tumors.

Mouse Monoclonal, Clone: Kp-1

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Cytoplasmic, Membranous



Tonsil

Catalog no.	Volume
C68-2010-1	0.1 ml concentrated
C68-2010-2	0.5 ml concentrated
C68-2010-3	1ml concentrated
C68-2010-4	1ml Prediluted
C68-2010-5	7 ml Prediluted
C68-2010	5 Positive control slides