




CD7

USA IVD

This is a useful marker in the identification of neoplastic proliferations.

Mouse Monoclonal, Clone: MRQ-12

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** IgG_{2b} **Cellular localization:** Membranous



Tonsil

Catalog no.	0.1 ml concentrated
CD7-2010-1	0.5 ml concentrated
CD7-2010-2	1 ml concentrated
CD7-2010-3	1 ml Prediluted
CD7-2010-4	7 ml Prediluted
CD7-2010-5	5 Positive control slides
CD7-2010	

Volume