




# CD74

USA IVD

This is used to differentiate atypical fibroxanthoma from malignant fibrous histiocytoma.

**Mouse Monoclonal, Clone: LN2**

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8<sup>0</sup>C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub> **Cellular localization:** Cytoplasmic, Membranous



Tonsil

<b>Catalog no.</b>	0.1 ml concentrated
C74-2010-1	0.5 ml concentrated
C74-2010-2	1ml concentrated
C74-2010-3	1ml Prediluted
C74-2010-4	7 ml Prediluted
C74-2010-5	5 Positive control slides
C74-2010	

## Volume