




CD79a

USA IVD

This is a B-cell marker that is often used to complement CD20.

Mouse Monoclonal, Clone: JCB117

Control: Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁/K **Cellular localization:** Membranous



Tonsil

Catalog no.	Volume
C79-2010-1	0.1 ml concentrated
C79-2010-2	0.5 ml concentrated
C79-2010-3	1ml concentrated
C30-2010-4	1ml Prediluted
C79-2010-5	7 ml Prediluted
C79-2010	5 Positive control slides