




## CD8

USA IVD

This is a useful marker for distinguishing helper-inducer T-lymphocytes, and most peripheral T-cell lymphomas.

**Mouse Monoclonal, Clone: C8/144B**

**Control:** Tonsil. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG<sub>1</sub>/K **Cellular localization:** Membranous



Tonsil

<b>Catalog no.</b>	0.1 ml concentrated
CD8-2010-1	0.5 ml concentrated
CD8-2010-2	1 ml concentrated
CD8-2010-3	1 ml Prediluted
CD8-2010-4	7 ml Prediluted
CD8-2010-5	5 Positive control slides
CD8-2010	

**Volume**