




## COX-2

USA IVD

This catalyzes the conversion of arachidonic acid to prostaglandin in H<sub>2</sub> in the first step in the biosynthesis of prostaglandin, thromboxanes, and prostacyclins.

**Rabbit Monoclonal, Clone: SP21**

**Control:** Colon Carcinoma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Cytoplasmic



Colon Carcinoma

<b>Catalog no.</b>	0.1 ml concentrated
Cox-2010-1	0.5 ml concentrated
Cox-2010-2	1ml concentrated
Cox-2010-3	1ml Prediluted
Cox-2010-4	7 ml Prediluted
Cox-2010-5	5 Positive control slides
Cox-2010	

**Volume**