




Caldesmon

USA IVD

This is used to label smooth muscle and tumors of smooth muscle, myofibroblastic, and myoepithelial differentiation.

Rabbit Monoclonal, Clone: E89

Control: Appendix . **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG **Cellular localization:** Cytoplasmic



Appendix

Catalog no.	Volume
CA-2010-1	
CA-2010-2	0.1 ml concentrated
CA-2010-3	0.5 ml concentrated
CA-2010-4	1 ml concentrated
CA-2010-5	1 ml Prediluted
CA-2010	7 ml Prediluted
	5 Positive control slides

