




Calponin-1

USA IVD

This can differentiate benign sclerosing lesions of the breast from carcinoma.

Rabbit Monoclonal, Clone: EP798Y

Control: Appendix. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8⁰C. **Shelf life:** Up to 36 months. **Isotype:** EP798Y, IgG, CALP, IgG1/K **Cellular localization:** Cytoplasmic



Appendix

Catalog no.	Volume
CB-2010-1	
CB-2010-2	0.1 ml concentrated
CB-2010-3	0.5 ml concentrated
CB-2010-4	1 ml concentrated
CB-2010-5	1 ml Prediluted
CB-2010	7 ml Prediluted
	5 Positive control slides