




# Calretinin

USA IVD

This can differentiate mesothelioma from adenocarcinomas of the lung and other sources.

## Rabbit Polyclonal

**Control:** Mesothelioma. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** **Cellular localization:** Cytoplasmic, Nuclear



Mesothelioma

Catalog no.	Volume
CD-2010-1	0.1 ml concentrated
CD-2010-2	0.5 ml concentrated
CD-2010-3	1 ml concentrated
CD-2010-4	1 ml Prediluted
CD-2010-5	7 ml Prediluted
CD-2010	5 Positive control slides

**Prediluted.** Smallest sizes available upon request