




Cytokeratin 19

USA IVD

This is used to identify thyroid carcinoma of the papillary type.

Mouse Monoclonal, Clone: A53-B/A2.26

Control: Colon. **Embedded:** Paraffin and frozen tissue. **Storage:** 2-8°C. **Shelf life:** Up to 36 months. **Isotype:** IgG₁ **Cellular localization:** Cytoplasmic



Colon

Catalog no.	0.1 ml concentrated
Y19-2010-1	0.5 ml concentrated
Y19-2010-2	1ml concentrated
Y19-2010-3	1ml Prediluted
Y19-2010-4	7 ml Prediluted
Y19-2010-5	5 Positive control slides
Y19-2010-	

Volume